



Technical Specifications

Product Name		D-Lux
Light Wavelength Range		420~490 nm
Light Intensity		1,600 mW/cm ²
Charging Time		Approx. 1 hrs 30 min
Operation time per full charge		160 times per 10 sec
Dimensions	Handpiece	165 x 165 x 25mm (L x W x H)
	Light Guide	11 x 11mm
	Light Guide Turbo (option)	11 x 8mm
	Charger	90 x 90 x 35mm (L x W x H)
Weight	Handpiece	125g
	Light Guide	30g
	Charger	150g
Power Supply (Medical Adaptor)	Input	100~240V AC/50~60 Hz
	Output	6V DC, 2.5A
Operating Conditions	Temperature	-10 ~ 40°C
	Atmospheric Pressure	800~1,060 (hpa)
	Humidity	10~85%
Transportation & Storage Conditions	Temperature	-10~45°C
	Atmospheric Pressure	600~1,060 (hpa)
	Humidity	10~90%

DiaDent Network

DiaDent Group International(Korea & International)
626, Yeonje-Ri, Gangoe-Myeon, Cheongwon-Gun, Chungcheong Buk Do, Rep. Of Korea
Tel(82-43-266-2315 Fax(82-43-262-9658
E-mail : diadent@diadent.co.kr

DiaDent Group International Inc.(Canada & USA)
11-3871 North Fraser Way, Burnaby B.C, Canada V5J 5G6
Tel(1-604-451-8851 Fax(1-604-451-8865
E-mail : diadent@diadent.com

DiaDent Europe B.V (Europe)
Antennestraat 70, 1322 AS Almere, The Netherlands
Tel(31-36-549-8607 Fax(31-36-536-7317
E-mail : diadent@diadent europe.com



©Copyright 2012 DiaDent Group Int
www.diadent.co.kr